

## INTRODUCTION

- o **bioSensAll™**: a powerful pharmacology BRET-based platform that enables the analysis and monitoring in real-time of signal transduction occurring in several pathways (GPCR or RTK (1)).
- o At Domain Therapeutics: development of a brand new and innovative pharmacological BRET platform dedicated to study Inhibitory Immune Checkpoint (ICP).
- o Inhibitory Immune Checkpoint (PD-1 or CTLA-4): front-runner targets in immunology field (2):
  - o involved in immune cells exhaustion,
  - o lead to a blockade of the anti-tumor response.
- o Blocking strategies targeting the PD-1/PD-L1 or CTLA-4 axis:
  - o proven and reliable therapies to strike cancer.
  - o more and more molecules such as blocking antibodies (Abs) or small molecule entities (SMEs) developed, and some of them are now in the drug market.


**In Drug Discovery, to screen and characterize these therapies: New robust and reliable cell-based assays are urgently needed.**

## OBJECTIVES AND METHODS


Development of BRET assays dedicated to either PD-1 or CTLA-4 axis:

- o Wild type forms of ICP receptors,
- o ICP receptors fused with rGFP,
- o co-culture with PD-L1/PD-L2 or CD80-CD86 positive cells to activate effector cells. Inspired from *in vivo* process where pathways are triggered by cancer cells expressing ICPs ligands.

### BRET Assays – PRINCIPLE




### BRET Assays - WORKFLOW



### OBJECTIVES


- PD-1 and CTLA-4 BRET assay validation based on:
- o monitoring the activation of PD-1 or CTLA-4 pathways in real time, and the specificity of the stimulation,
  - o assessing the inhibitory activity of know ICPs inhibitors (SMEs or blocking antibodies.)

## Activation of PD-1 or CTLA-4 pathways in BRET assays




**Figure 1. PD-1 and CTLA-4 BRET assays respond to ligand stimulation.**

For both targets and both assays, a BRET ratio is observed in presence of ligand expressing cells. This activation of the pathway is dependent of the ligand quantity.



**Figure 2. PD-1 BRET assays respond to PD-L2 stimulation.**


Upon binding to its ligands PD-L1 or PD-L2, PD-1 or PD-1-rGFP engaged SHP2.



**Figure 3. The phosphodeficient version of the key amino acid for SHP2 recruitment inhibits BRET signal.**

As previously characterized (3), Y248F point mutation of PD-1 is not able to recruit SHP2 phosphatase. This is perfectly translated in our PD-1-rGFP assay by a complete loss of the BRET signal.


## Monitoring of pathways activation in real time



**Figure 4. Kinetic of the pathway activation in PD-1 and CTLA-4 assays.**

The kinetic activation of PD-1 and CTLA-4 pathways can be monitored in real-time. The maximal BRET signal is reached at 1H of incubation.


## Specificity of pathway activation - PD-1



**Figure 5. Only PD-L1 cells can stimulate PD-1 BRET assays.**

CD86 positive cells are not able to stimulate PD-1 BRET assay.


## Specificity of pathway activation - CTLA-4



**Figure 6. Only CD80-86 cells lead to an activation of the CTLA-4 BRET assays.**

PD-L1 positive cells are not able to stimulate CTLA-4 assays.


## PD-1 BRET assays validation with SMEs




**Figure 7. Pharmacology profile of PD-1/PD-L1 inhibitors into PD-1 assays.**

Both assays are able to highlight the inhibitory effect of SMEs described to block PD-1 pathway: BMS-1166 and BMS-1001. IC50 are similar compared to (4).


Interestingly, the assays can reveal the specificity of the ligand binding as no inhibitory effect is observed with a co-culture of PD-L2 expressing cells.



### Bosutinib, inhibitor of Src kinase family




### CA-170



No inhibition of BRET signal was observed in the 2 assays. As already characterized (3,4), this molecule is not able to interact with PD-1/PD-L1 interface.


## PD-1 BRET assays validation using PD-1 blocking antibody



**Figure 8. Pharmacology profile of well characterized PD-1 antibody: Nivolumab into PD-1 assays.**

PD-1 neutralizing antibody Nivolumab inhibits BRET signal of PD-1 assays with PD-L1 or PD-L2 stimulation. BRET assays highlight the blocking activity of this antibody in both stimulations (PD-L1 and PD-L2).

## CTLA-4 BRET assays validation using blocking CTLA-4 antibody



**Figure 9. Pharmacology profile of CTLA-4 neutralizing antibody into CTLA-4 assays.**

The two assays highlight the inhibitory effect of Ipilimumab antibody. As a specific control, Nivolumab doesn't inhibit CTLA-4 pathway.

## CONCLUSION

**bioSensAll™** is a powerful pharmacological platform:

- o A versatile pharmacological platform able to be highly adaptable for hot topic targets such as **Immune Checkpoint targets**:
  - o **Primary and secondary assays** for PD-1 and CTLA-4 axis were developed,
  - o Assays were able to reveal the inhibitory activities of therapeutic molecules (SMEs, mAbs)
- o **Adaptably of the platform across several biosensors,**
- o Capable to monitor signaling pathways in real time kinetics for inhibitory ICPs, and will be assessed for stimulatory ICPs,
- o Can be companion tests to enhance drug discovery and also can be used as **quality control/batch release testing** for therapeutic molecules.

## REFERENCES

1. <http://biosensall.com/biosensall/>
2. Axel Hoos et al., *Nature Reviews Drug Discovery*, 2016
3. Michael Peled et al., *PNAS*, 2017
4. Lukasz Skalniak et al., *Oncotarget*, 2017
5. Aravindhan Ganesan et al., *Scientific Reports*, 2019
6. Bogdan Musielak et al., *Molecules*, 2019, ver